

Clench[®] Compression Staple



Sterile - Single Use Nitinol Compression Staple System

7x5



9x7



11x10



15x12



18x18




20x20




Technical Features

- No heat activation required
- Constant compression
- Smooth profile
- All-in-one instruments & implant kit
- Sterile, single use kits
- Reloadable instrument design
- Multiple tamp options
- Robust instrumentation



SIZE	WIDTH (W)	LENGTH (L)	DIAMETER	COMPRESSION
7mm x 5mm	7mm	5mm	1.65mm	0.50mm
9mm x 7mm	9mm	7mm		0.63mm
11mm x 10mm	11mm	10mm		1.15mm
13mm x 10mm	13mm	10mm		1.15mm
15mm x 12mm	15mm	12mm	2.25mm	1.25mm
15mm x 15mm	15mm	15mm		1.60mm
18mm x 15mm	18mm	15mm		1.60mm
18mm x 18mm	18mm	18mm		2.15mm
20mm x 20mm	20mm	20mm		2.25mm










PART #	CONTENTS
152-40001	(1) Implant Sizer

PART #	SIZE	CONTENTS
152-30705	7mm x 5mm	(1) Drill Guide (1) Drill
152-30907	9mm x 7mm	(2) Guide Post (1) Delivery Instrument (1) Nitinol Staple
152-31110	11mm x 10mm	
152-31512	15mm x 12mm	Note: Contents are specific to the implant size included. Do not use instruments with other implant sizes.
152-31818	18mm x 18mm	
152-32020	20mm x 20mm	

PART #	SIZE	CONTENTS
152-40705	7mm x 5mm	(1) Delivery Instrument (1) Nitinol Staple
152-40907	9mm x 7mm	
152-41110	11mm x 10mm	
152-41512	15mm x 12mm	
152-41818	18mm x 18mm	
152-42020	20mm x 20mm	

Surgical Procedure



This brochure contains information intended to highlight a Reign Medical product. It is highly recommended to review the package insert, surgical technique guide, and/or participate in a product training course before using any Reign Medical product. Surgeons should use their medical experience when choosing which products and techniques to use for each patient. Please contact your authorized representative for any questions and product availability.



6303 E. 102nd Street
Tulsa, OK 74137
www.reignmedical.com

Office: 918-302-9546
Fax: 918-844-4356
© 2018 Reign Medical, LLC. All Rights Reserved.